FICHE TECHNIQUE

Hemostatic gauze - Gaze hémostatique - Soluble - 1-Medical

Designed to quickly absorb exudates from the wound surface, forming a viscous gel that seals the wound for trauma and post-operative wounds.

The Soluble Hemostatic Gauze is designed to quickly absorb exudates from the wound surface, seal the wound, and accelerate healing for trauma and post-operative wounds. When the gauze absorbs exudates, creating a gelatinous mass that fills wound voids, seals capillary ends, and aids in wound healing.

Technical specifications :

- 35,4 x 3,9 in - 900 x 100 mm,
- Weight : 20 g,
- Primarily made of sodium carboxymethyl cellulose,

Instruction for use :

- Remove it from its sterile packaging and cut it to fit the wound area,
- Apply slight direct pressure until the bleeding stops (the dosage can be adjusted based on the wound needed), moistening with saline for better adherence,
- If desired, combine with standard wound dressings for additional protection.
- If the gauze becomes saturated, moisten it with sterile water or saline solution, then gently remove.

Indications and Contraindications :

The Soluble Hemostatic Gauze is ideal for external, short-term use in controlling bleeding in a general population, including those on anticoagulant medications. However, it is not intended for internal use.

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