

## FICHE TECHNIQUE

### Assault medical backpack - Sac à dos médical d'assaut - R20 - 1-Medical

Complete backpack to support tactical teams during their missions on the field

This backpack is a complete and robust backpack to empower tactical teams during their missions. Perfectly equipped to carry a large volume of emergency equipment and supplies for advanced cardiac and trauma life support, this backpack is your trusty companion during deployment or critical incidents.

#### Technical specifications :

- Dimensions : 18 x 12 x 5 in - 460 x 300 x 130 mm,
- Empty weight : 1,7 kg,
- Capacity : 20 L,
- Material : Nylon 1000 D.

#### Configuration :

3x chest seal twin pack,  
5x compressed gauzes,  
3x 4" emergency bandage,  
3x 6" emergency bandage,  
5x tourniquet,  
2x 36" aluminium splint,  
1x medical shears,  
3x nasal airway,  
3x lubricant jelly,  
2x soluble hemostatic gauze,  
1x multifunctional emergency shears,  
1x reusable led penlight,  
10x gloves,  
3x combat cravat,  
1x resuscitator,  
5x marker,  
6x alcohol pad,  
5x CPR face shield,  
5x combat tape,  
5x emergency blanket,  
6x eye shield,  
5x elastic bandage,  
4x oroharyngeal airway,  
5x burn dressing (3,5g),  
2x burn dressing (100\*100mm),  
1x burn dressing (200\*450mm).

CE medical equipment.

Ce sac à dos est un sac complet et robuste conçu pour soutenir les équipes tactiques lors de leurs missions. Parfaitement équipé pour transporter un grand volume d'équipement d'urgence et de fournitures pour les soins avancés en traumatologie et en réanimation cardiaque, ce sac à dos est votre fidèle compagnon lors du déploiement ou des incidents critiques.

#### Spécifications techniques :

Dimensions : 18 x 12 x 5 pouces - 460 x 300 x 130 mm

Poids à vide : 1,7 kg



Capacité : 20 L  
Matériau : Nylon 1000 D

Configuration :

3x paquets jumeaux de joints thoraciques  
5x gazes compressées  
3x bandages d'urgence 4"  
3x bandages d'urgence 6"  
5x garrots  
2x attelles en aluminium 36"  
1x ciseaux médicaux  
3x canules nasales  
3x gel lubrifiant  
2x gazes hémostatiques solubles  
1x ciseaux d'urgence multifonctionnel  
1x lampe de poche à LED réutilisable  
10x gants  
3x cravates de combat  
1x reanimateur  
5x marqueurs  
6x tampons d'alcool  
5x masques de réanimation cardio-pulmonaire  
5x ruban de combat  
5x couverture de survie  
6x écrans oculaires  
5x bande élastique  
4x canules oro-pharyngées  
5x pansement pour brûlures (3,5 g)  
2x pansement pour brûlures (100\*100 mm)  
1x pansement pour brûlures (200\*450 mm)

Matériel médical CE.