

**Compressed gauze - Gaze comprimée -
Quick-response - 1-Medical**

6-ply absorbent cotton gauze for the treatment of bleeding wound.

This compressed gauze is a sterile 6-ply woven absorbent cotton gauze for the treatment of bleeding wounds. It can be pulled out from the top, to directly fill the wound.

Technical specifications :

- Size : 4,5 in x 4,1 yds - 114 mm x 3,75 m,
- Weight : 41 g,
- Material : 100% cotton.

CE medical equipment.



Cette gaze compressée est une gaze de coton absorbante stérile à 6 plis pour le traitement des plaies saignantes. Elle peut être déroulée par le haut pour remplir directement la plaie.

Spécifications techniques :

- Taille : 4,5 pouces x 4,1 yds - 114 mm x 3,75 m,
- Poids : 41 g,
- Matériau : 100 % coton.

Matériel médical CE.