

### Decompression needle - Aiguille de décompression - 1-Medical

The decompression needle is an emergency air release to manage pneumothorax.

The decompression needle is a vital and indispensable tool designed for emergency situations. When every second counts, the decompression needle equips medical professionals, emergency responders, and military medics to manage tension pneumothorax injuries in the field effectively.

#### Technical specifications :

- Decompression needle - ADV-14Ga : 14 ga x 3,25 in (82 mm),
- Decompression needle - ADV-10Ga : 10 ga x 3,25 in (95 mm),
- Weight : 18 g.

CE medical equipment.



---

L'aiguille de décompression est un outil vital et indispensable conçu pour les situations d'urgence. Lorsque chaque seconde compte, l'aiguille de décompression équipe les professionnels de la santé, les secouristes d'urgence et les médecins militaires pour gérer efficacement les blessures de pneumothorax sous tension sur le terrain.

#### Spécifications techniques :

- Aiguille de décompression - ADV-14Ga : 14 G x 3,25 pouces (82 mm),
- Aiguille de décompression - ADV-10Ga : 10 G x 3,25 pouces (95 mm),
- Poids : 18 g.

*Matériel médical CE.*