

Emergency thermal blanket - Couverture de survie - 1-Medical

Ideal for emergencies that require warmth or protection from all-weather conditions.

These Emergency Thermal Blankets are an essential piece of any emergency kit, perfect for bug-out bags, go-bags, and home/auto emergency kits.

Constructed from high-quality, waterproof, and weatherproof foil and Mylar, these blankets ensure heat retention or deflection. These blankets are a lifeline in damp, wet, and cold environments, where they block rain, snow, and moisture.

Technical specifications :

- Size unfolded : 63 x 82,7 in (160 x 210 cm),
- Weight : 40 g.

CE medical equipment.

Ces couvertures thermiques d'urgence sont un élément essentiel de tout kit d'urgence, parfaites pour les sacs de secours, les sacs de voyage d'urgence et les kits d'urgence à domicile ou pour la voiture.

Fabriquées à partir de feuille de qualité supérieure imperméable et résistante aux intempéries, ces couvertures garantissent la rétention ou la réflexion de la chaleur. Elles sont une bouée de sauvetage dans les environnements humides, humides et froids, où elles bloquent la pluie, la neige et l'humidité.

Spécifications techniques :

- Taille dépliée : 63 x 82,7 pouces (160 x 210 cm),
- Poids : 40 g.

Matériel médical CE.

