

Multi-purpose emergency litter - Brancard pliable - 1-Medical

Emergency foldable and space-saving litter for emergency situation.

An emergency carrying solution designed for swift deployment in mass casualty incidents, natural disasters, active shooter situations, and routine patient transport.

Technical specifications :

- Stretcher size : 73 x 27 in - 1860 x 680 mm,
- Package size : 8,3 x 12,6 in - 210 x 320 mm,
- Weight : 582 g,
- Load rating : 206 kg,
- Material : Nylon 500d, Water and tear-resistant material
- Ten large, robust carrying handle.

Space-saving foldable design for easy storage and portability, includes a storage bag for convenient protection and transport of the product.

CE medical equipment.

Une solution de transport d'urgence conçue pour un déploiement rapide lors d'incidents de masse, de catastrophes naturelles, de situations d'attaques actives et de transport de patients de routine.

Spécifications techniques :

- Taille du brancard : 73 x 27 pouces - 1860 x 680 mm,
- Taille de l'emballage : 8,3 x 12,6 pouces - 210 x 320 mm,
- Poids : 582 g,
- Capacité de charge : 206 kg,
- Matériau : Nylon 500d, matériau résistant à l'eau et aux déchirures
- Dix grandes poignées de transport robustes.
- Design pliable peu encombrant pour un rangement et un transport faciles, comprend un sac de rangement pour une protection et un transport pratiques du produit.

Matériel médical CE.

