



### Tactical IFAK - Kit tactique de premiers secours - 002 Mini - 1-Medical

Versatile and adaptable IFAK to cater to multiple needs

This IFAK is a versatile and adaptable first aid kit designed to cater to the unique needs of various users. Whether you're a first responder, outdoor enthusiast, or a professional, this customizable IFAK is there to support you during emergency situations.

#### Technical specifications :

- Dimensions : 6,3 x 3,5 x 2,8 in - 160 x 90 x 70 mm,
- Empty weight : 100 g,
- Weight with configuration : 113 g,
- Material : Nylon 500 D.

#### Configuration :

- 1x 4" emergency bandage,
- 1x gloves,
- 1x scissors,
- 1x tourniquet,
- 1x compressed gauze.

CE medical equipment.

---

Cet IFAK est un kit de premiers secours polyvalent et adaptable conçu pour répondre aux besoins uniques de différents utilisateurs. Que vous soyez un premier intervenant, un passionné de plein air ou un professionnel, cet IFAK personnalisable est là pour vous soutenir lors de situations d'urgence.

#### Spécifications techniques :

- Dimensions : 6,3 x 3,5 x 2,8 pouces - 160 x 90 x 70 mm
- Poids à vide : 100 g
- Poids avec configuration : 113 g
- Matériau : Nylon 500 D

#### Configuration :

- 1x bandage d'urgence de 4 pouces
- 1x paire de gants
- 1x ciseaux
- 1x garrot
- 1x gaze compressée

Matériel médical CE.