

Trauma scissors - Ciseaux d'urgence - 1-Medical

High-Quality, Multi-Functional Shears

Scissors with sharp, 18.5 cm long stainless steel blades featuring milled serrations, these shears can cut through various hard objects with ease.

The design is topped off with injection-molded handles and a steel-reinforced carabiner, making these shears not just highly functional, but durable and dependable.

Technical specifications :

- Size : 18,5 cm,
- Color : black.

CE medical equipment.



Ces ciseaux sont équipés de lames en acier inoxydable tranchantes de 18,5 cm de long avec des dentelures usinées, leur permettant de couper facilement à travers divers objets durs.

Le design est complété par des poignées moulées par injection et un mousqueton renforcé en acier, ce qui rend ces ciseaux non seulement très fonctionnels, mais aussi durables et fiables.

Spécifications techniques :

- Taille : 18,5 cm,
- Couleur : noir.

Matériel médical CE.