

FICHE TECHNIQUE

Triangular bandage - Bandage triangulaire - 1-Medical

Large and easy to tear triangular bandage

Large triangular bandage to use on the field. For more praticicity this bandage is easy to tear which allow a wide variety of uses.

Technical specifications :

- Dimensions : 38 x 38 x 53 in - 965 x 965 x 1346 mm,
- Weight : 37 g,
- Material : cotton.

CE medical equipment.



Grand bandage triangulaire à utiliser sur le terrain. Pour plus de praticité, ce bandage est facile à déchirer, ce qui permet une grande variété d'utilisations.

Spécifications techniques :

- Dimensions : 38 x 38 x 53 pouces - 965 x 965 x 1346 mm,
- Poids : 37 g,
- Matériau : coton.

Matériel médical CE.