

**Vented chest seal - Joint thoracique ventilé -
Advanced - 1-Medical**

The advanced emergency solution to prevent pneumothorax in case of open chest wound.

Chest seals are the emergency solution for open chest wounds. They are designed to prevent pneumothorax, the second most common cause of preventable death on the battlefield.

Technical specifications :

- Dimensions : 4 x 6 in/pc - 105 x 150 mm/pc,
- Weight : 40 g/pc,
- Material : hydrocolloid,
- Improve sealing performance,
- Reduce fluid reflux,
- Multi layers

CE medical equipment.

Les pansements thoraciques sont la solution d'urgence pour les plaies ouvertes à la poitrine. Ils sont conçus pour prévenir le pneumothorax, la deuxième cause de décès évitable la plus courante sur le champ de bataille.

Spécifications techniques :

- Dimensions : 4 x 6 pouces/pc - 105 x 150 mm/pc,
- Poids : 40 g/pc,
- Matériau : hydrocolloïde,
- Améliore les performances d'étanchéité,
- Réduit le reflux de fluide,
- Multi-couches.

Matériel médical CE.

